

"We respect Dentistry, Dental Technology and Science which makes us stand out from our competition"
Kash Qureshi, Managing Director

All Ceramic ... Zirconia VS E.max, Which one?

All Ceramics table - IPS e.max VS Zirconia? Clinical indications, flexural strength. A simple table to for treatment planing & giving treatment options chair side.

All Ceramic Material	Flexural Strength	Clinical Indications	Aesthetics or Function
Lithium Disilicate IPS e.max	360 MPa to 400 MPa	Inlays, onlays, 3/4 crowns, single unit crowns, short span FPD's in anterior and premolar region, implant abutments.	Aesthetics
Zirconia Cercon	900 MPa - 1,200 MPa	Single unit anterior crowns, single unit posterior crowns, Large span anterior and posterior bridges, implant abutments, crowns & bridges.	Function
Alumina In-ceram Alumina Procera	450 Mpa 600 MPa	Single unit crowns. Single unit crowns, anterior bridges.	Aesthetics
Leucite IPS empress	Leucite IPS empress	Veneers, inlays, onlays, 3/4 crowns, single unit crowns.	Aesthetics & Function




Full Zirconia



Cercon (Zirconia)



Smile Make Overs



IPS e.max

**Call or Email
today to discuss
a case that may
be suitable:**

0208 520 8520

office@bremadent.co.uk

www.bremadent.co.uk

www.facebook.com/bremadent

0208 520 8528
www.bremadent.co.uk
office@bremadent.co.uk




"We respect Dentistry, Dental Technology and Science which makes us stand out from our competition"
Kash Qureshi, Managing Director

APPRENTICES WANTED?

Do you know any school leavers that are looking for a career in Dentistry? Contact Kash for further information on our apprentice and training scheme for young persons.

We often get asked what does Bremadent look like, here is a gallery of our Dental Laboratory...



Are you ready for the future of Dentistry?



The TRIOS 3D intraoral scanner enables you to create accurate digitally impressions and validate them instantly chair-side.

Benefits for the Dentist:

Consistent

- High quality accuracy from digital scans
- Instant preparation and impression validation
- Improved accuracy and clinical results
- Better fitting restorations
- Send to Bremadent instantly online

Productive

- Faster than traditional impression taking
- Less adjustment and grinding during seating
- No impression materials or mess
- Reduced need for retakes
- Comfortable & Quick

Cost effective

- Reduced over-all chair-time
- Cost savings related to materials
- Improved clinical results



CONTACT BREMADENT FOR MORE INFO!