

0208 520 8528
www.bremadent.co.uk
office@bremadent.co.uk



“We respect Dentistry, Dental Technology and Science which makes us stand out from our competition”
Kash Qureshi, Managing Director

Start of a New Year!

Bremadent was established in 1969 as a small one level, crown & bridge laboratory adjacent to a Dental Practice. In 2017, Bremadent is a full service Dental Laboratory, two storeys high with over 10 departments and in-house fully functional denture clinic for Patient's and Dentist's.

We offer Implantology, Crown & Bridge, Prosthetics, Orthodontics, Vacuum formed appliances and a mixture of combination work for Private & NHS Dentist in London and the UK with the help of the digitalisation of Dental Technology.

Why not kick start your year with a Dental Laboratory with history.....




How much should I charge for Dentures?

Below is a helpful chart to show the average fee for denture patients in accordance with location, how does your fee compare with the national average?

<http://www.dentistry.co.uk/2016/11/24/dentistry-fee-s-survey-results-dentures/>

Fun Dental Fact!

Everyone's set of teeth are unique to them, just as fingerprints are. Even identical twins can be told apart from their teeth.



0208 520 8528
www.bremadent.co.uk
office@bremadent.co.uk



"We respect Dentistry, Dental Technology and Science which makes us stand out from our competition"
Kash Qureshi, Managing Director



Here is a list of key components that make up a Smile Make Over case:



- Study Models
- Denar Articulation
- Preparation guide models
- Diagnostic wax ups
- Vacuum Formed Stent
- Shade taking
- Digital Photography
- IPS e.max or Zirconia restorations

Whether you're on a course or have a complex case, Bremadent can help. Call us on 0208 520 8528.

Clinical Indications for IPS e.max

- Inlays, onlays
- Veneers
- 3/4 crowns
- single unit crowns
- short span FPD's in anterior and premolar region
- implant abutments

IPS e.max the Lithium Disilicate glass ceramic with a flexural strength of over 400mpa. Speak to our Ceramics department to find out more!



Download our **NEW!** app today!

SEARCH: BREMADENT



What is a Kois Appliance?

A Kois deprogrammer appliance is designed for patients with TMJ disorders, bruxism, headaches & pain, muscle fatigue, occlusal dysfunction, protection against wear and establishing centric relation as a diagnostic tool. When mastication is performed the TMJ adjusts automatically to accommodate and creates a pattern.

De-programmers work when anterior teeth are in contact and the muscle contraction is forced. This causes the lateral pterygoid muscles to relax and allows the TMJ to seat in its most comfortable and natural position. It allows muscle relaxation as it removes occlusal interferences, habitual functional movements the patients has developed to avoid pain.

A Kois de-programmer is a removable device that covers the palatal aspect of the maxilla with an anterior incisal platform between the upper central's palatally. The platform follows the occlusal plane of the arch and separates the maxilla and mandible this creates a contact between the platform and the lower central's. Protocols recommend a labial arch wire or two clasps posteriorly are added for retention. The device can be worn for up to 20 hours a day (except meals) or at night and is worn until the necessary muscle deprogramming is accomplished which can be days or weeks.
Contact Kash today to find out more: 0208 520 8528